

----- PROCESSING PARAMETERS -----  
dc\_balance( 0, FALSE )  
seXp( 0.2[Hz], 0.0[s] )  
trapezoid( 0[%], 0[%], 80[%], 100[%] )  
zeroFill( 1 )  
fft( 1, TRUE, TRUE )  
machinephase  
ppm  
auto\_reference( 5[%], TRUE )  
auto\_reference( 5[%], TRUE )

Filename = Ace 8\_non\_1H-1-7.jdf  
Author = delta  
Experiment = single\_pulse.jxp  
Sample\_Id = Ace 8  
Solvent = CHLOROFORM-D  
Creation\_Time = 6-OCT-2015 15:11:43  
Revision\_Time = 6-OCT-2015 14:33:59  
Current\_Time = 19-OCT-2015 11:06:16

Data\_Format = 1D COMPLEX  
Dim\_Size = 13107  
Dim\_Title = Proton  
Dim\_Units = [ppm]  
Dimensions = X  
Site = ECS 400  
Spectrometer = DELTA2\_NMR

Field\_Strength = 9.389766[T] (400[MHz])  
X\_Acq\_Duration = 2.18365952[s]  
X\_Domain = 1H  
X\_Freq = 399.78219838[MHz]  
X\_Offset = 5[ppm]  
X\_Points = 16384  
X\_Prescans = 1  
X\_Resolution = 0.45794685[Hz]  
X\_Sweep = 7.5030012[kHz]  
X\_Sweep\_Clippped = 6.00240096[kHz]  
Irr\_Domain = Proton  
Irr\_Freq = 399.78219838[MHz]  
Irr\_Offset = 5[ppm]  
Tri\_Domain = Proton  
Tri\_Freq = 399.78219838[MHz]  
Tri\_Offset = 5[ppm]  
Clipped = FALSE  
Scans = 8  
Total\_Scans = 8

Relaxation\_Delay = 5[s]  
Recvr\_Gain = 58  
Temp\_Get = 23.9[degC]  
X\_90\_Width = 12[us]  
X\_Acq\_Time = 2.18365952[s]  
X\_Angle = 45[deg]  
X\_Atn = 2.4[dB]  
X\_Pulse = 6[us]  
Irr\_Mode = Off  
Tri\_Mode = Off  
Dante\_Presat = FALSE  
Initial\_Wait = 1[s]  
Repetition\_Time = 7.18365952[s]

